

Optimizing patient outcomes in Cataract and Refractive Lens Exchange Surgery



The Company

EyePCR is a medical technology company focused on eliminating the post-operative complications experienced by 30-40% of cataract surgery patients

Patented fixOflex is a unique intraocular ophthalmic platform designed to preserve capsular form and intracapsular space similar to its natural form and address complications associated with cataract and refractive lens exchange surgery

CE Mark certification

fixoflex received CE Mark certification under the European Union Medical Device Regulation (EU MDR 2017/745) In early 2026



Main Indications of use:

- Preservation of open capsule and intracapsular space
- PCO prevention (up to at least 12 months post op)

Peer reviewed publications

Pallikaris et al, TVST, 2026;15(2):8

- AEs within the range of published data and ISO, all resolved
- Reduced number of PCO cases by 93%
- No need for YAG laser capsulotomy in any of the patients (vs. 3 in control group)

Britz et al, TVST, Accepted and pending publication

- Optical Quality Preservation
- Stable fixation of IOL designs within fixOflex

EYEPCR is launching fixOflexcommercially



Japan
Launch 2024 (in market already)

Europe
(Germany, Spain, Italy)
Singapore

Q2-3 2026



End of Q4 2026

European Markets
Australia
Egypt

India
Latin America
Asia

2027-2028



Post-FDA

USA - FDA
Regulatory Plan
in Progress

Raising € 4.0 MM for EU, Asia and EM expansion

Addressing the largest number of global surgeries

Unique solution ready for a growing market

30+
Million

Cataract
Surgeries
Worldwide


\$2.5+
Billion

Potential
Global
Market

3x

Faster
growth for
pIOL
surgeries





Alexandria Main Clinical Study : 3 Years Results, 68 patients



■ fixOflexTM + Hydrophobic IOL ■ Hydrophilic IOL

- 77% of patients with just a Hydrophilic IOL had PCO, whereas only **4.5% (not clinically significant, Grade 1-3)** of fixOflex patients were observed with mild PCO.
- **18.8%** of patients from the control group (**Grade 4 PCO**) undergone Nd:YAG capsulotomy in comparison to 0% of the fixOflex patients.

Japanese Study with premium IOL: 3 Months Results, 38 patients per group



Bilateral implantation combined with a trifocal IOL, Control group


- Mean Binocular UDVA and UNVA statistically better in the device group with a higher percentage of patients in the device group achieving binocular 20/25 or better
- Significantly higher CS under mesopic or glare conditions in the device group
- VFQ-25 scores significantly higher scores in the peripheral vision subdomain in the device group

Athens Pediatric Study (Anticipated 30 patients) : Case Report, 1 Year follow up



Contralateral Comparative Study (fixOflex+IOL vs Conventional IOL Implantation)

- 13-year-old boy with bilateral congenital cataract, first diagnosed at the age of 10 years.
- Preoperative visual acuity: 20/100 in the right eye and 20/80 in the left eye.
- Right eye: standard cataract surgery with intraocular lens implantation.
- Left eye: implantation of the fixOflex ring prior to intraocular lens insertion.




RF: +1.00 -1.00 x135

UDVA: 20/32

BCVA: 20/20

Slit-lamp: clear visual axis, without signs of PCO or anterior capsule contraction.



RF: +2.00 -1.00 x5

UDVA: 20/50

BCVA: 20/32



Slit-lamp: PCO grade 2 (EPCO) and significant anterior capsule contraction.

Benchmark Studies: Dysphotopsia and Stress



Dysphotopsia

Determinants of Zonular Stress During Intraocular Implantation

Courtesy of L. Britz, M. Hammer,, HS Son, and GU Auffarth
University of Heidelberg, Submitted for publication




fixoflex suppressed the characteristic crescent associated with negative dysphotopsia



Peak capsular bag expansion during implantation for each device: (A) fixOflex, (B) C-loop IOL, and (C) plate-haptic IOL

Capsular bag expansion induced by the leading haptic along the injection vector (white arrows).

While C-loop IOL exhibits a predominant expansion axis, capsular expansion during fixOflex implantation is distributed in a semicircular pattern across the capsular bag.



Capsular bag occupancy areas for (A) fixOflex, (B) plate-haptic IOL and (C) one-piece C-loop IOL.

Long axis (white, L) - Short axis (black, S)

fixOflex device demonstrated a consistently high and near-total capsular bag occupancy
The fixOflex showed a long/short ratio of 0.99 (0.97–1.00), the plate-haptic a ratio of 1.04 (1.00–1.14), and the one-piece C-loop of 1.10 (1.09–1.20).